

Our highest resolution 3D imaging, faster.

3DQuorum™ Imaging Technology, powered by Genius AI™

3DQuorum™ technology utilises Genius AI™-powered analytics to uniquely reconstruct high-resolution 3D Mammography™ data to produce 6mm SmartSlices. SmartSlices are designed to expedite reading time by reducing the number of images to review, with no compromise in image quality, sensitivity or accuracy.¹



Genius AI™-powered software designed to provide an accurate and timely diagnosis

Genius AI™ identifies clinically relevant regions of interest and preserves important features during the creation of SmartSlices. Each high-resolution SmartSlice overlaps the previous one by 3mm, ensuring there is no loss of 3D™ image data and continuity in scrolling.¹




Designed to expedite Hologic Clarity™ HD high-resolution 3D Mammography™ exam read times

When a Radiologist reads SmartSlices instead of 1mm slices, the number of 3D™ images to review is reduced by 2/3, leading to an average interpretation time saving of 1 hour per day, based on 8 hours of image interpretation time per day.¹



High-Resolution file size reduction

3DQuorum™ technology reduces the typical Hologic Clarity HD and Intelligent 2D™ study size by over 50%, bringing the storage space and network impact back down to that of standard resolution 3D™ imaging.



Images shown are high-resolution and for display purposes only.

Product Information

Minimum Hardware/Software Requirements

To ensure optimal workflow of 3DQuorum™ SmartSlices, the minimum requirements for Dimensions systems and connected Hologic products are listed below. For customers using PACS, it is recommended that customers work with their individual PACS vendor to ensure a successful integration.

	Selenia® Dimensions® / 3Dimensions™	Advanced Workflow Manager (AWM)	Reading Workstation	CAD & Quantra™	SecurXchange®
Software	Selenia® Dimensions®: <ul style="list-style-type: none"> Software v1.10 Hologic Clarity HD 3D™ Imaging 3Dimensions™: <ul style="list-style-type: none"> Software v2.1 Hologic Clarity HD 3D™ Imaging 	Software v1.10	SecurView®: v10.4 Unifi™ Workspace: v1.0.	Cenova™: 3.0 ImageChecker® CAD: 10.0 Quantra™: v2.2	Latest Available
Computer / Hardware	ASY-13349 <ul style="list-style-type: none"> For both Selenia® Dimensions® and 3Dimensions™ 	CMP-01377 and CMP-01378	Workstation: <ul style="list-style-type: none"> Dell Precision T3620 Video Card: <ul style="list-style-type: none"> MXRT-7600 Monitor: All currently serviced <ul style="list-style-type: none"> 12MP Coronis Uniti 5MP Tomo (5221, 5121, 6121, 5121-MA) SecurView® Manager: Dell PowerEdge T630 or R730	CMP-00712 and above	Dell PowerEdge R730
Operating System	Windows 10	Windows Server 2016	Workstation: Windows 10 SecurView® Manager: Windows Server 2016	Windows 7	LINUX

Ordering details

Part Number	Description
DIM-LIC-3DQ	3DQuorum™ license for new Selenia® Dimensions® or 3Dimensions™ system
DIM-LIC-3DQ-UP	3DQuorum™ license upgrade option for existing Selenia® Dimensions® or 3Dimensions™ system
DIM-LIC-TRIAL-3DQ	3DQuorum™ 90 days trial license for Selenia® Dimensions® and 3Dimensions™ systems with Clarity HD™

References

*Compared to 2D alone

1. Data on File: Clinical Study Report CSR-00116

SSL-00369-CAN-EN Rev 001 (01/22) Hologic Inc. ©2022 All rights reserved. Hologic, 3D, 3D Mammography, 3Dimensions, 3DQuorum, Cenova, Dimensions, Genius AI, Hologic Clarity, ImageChecker, Intelligent 2D, Quantra, SecurView, SecurXchange, Selenia, Unifi and The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

www.3dimensionsmammography.eu | canada2@hologic.com | 1-877-209-7192

No Compromise. No Comparison.