

# Reveal more with low-dose, highly accurate 3D Mammography™ exams.

## C-View™ Software

Raise your breast cancer screening performance<sup>1-7</sup> while minimizing patient radiation and discomfort with C-View™ software. An advanced algorithm takes high-quality tomosynthesis data and instantly generates 2D images to enhance details and speed analysis.



**Deliver superior clinical performance than 2D mammography alone for all breast types.<sup>1-7</sup>**



**Reduce the risk of retakes with ultra-fast 3.7-second scans.<sup>1-7</sup>**



**Greater accuracy at a lower dose.<sup>1, 4-6, 8-9</sup>**


**C-View™**  
Software



# The proof is in the details.

C-View 2D images are clinically proven and FDA approved to diagnostically replace the FFDM images within a tomosynthesis screening exam. The images also serve as a navigational aid to the tomosynthesis slice review. Published studies show that the low dose 3D Mammography™ exam finds invasive cancer earlier,<sup>1-4</sup> compared to 2D alone, while also reducing false positive recall rates.<sup>1,4-6</sup>

**See more, do more.**



**Architectural distortions, mass lesions, and bright spots commonly found in microcalcifications are more visible in the C-View 2D image than on the traditional FFDM 2D image or tomo slice.**

## Product information

C-View 2D imaging is a purchasable option available on both Selenia® Dimensions® and 3Dimensions™ systems. It is compatible with standard resolution 3D™ imaging only- not compatible with high-resolution 3D™ imaging. Refer to Dimensions Platform datasheet for additional technical product information.

## Imaging Modes

Combo Mode	Standard resolution 3D™ imaging + FFDM
TomoHD Mode	Standard resolution 3D™ imaging + C-View
ComboHD Mode	Standard resolution 3D™ imaging + FFDM + C-View

## Ordering details

Part Number	Description
SDM-LIC-0005	C-View 2D imaging software license

## References

**1** FDA PMA submission P080003/S001 physician labeling **2** Skaane P, Bandos AI, Eben EB, et al. Two-view digital breast tomosynthesis screening with synthetically reconstructed projection images: comparison with digital breast tomosynthesis with full-field digital mammographic images. *Radiology*. 2014 Jun;271(3):655-63. **3** Zuley M, Guo B, Catullo V, et al. "Comparison of Two-dimensional Synthesized Mammograms versus Original Digital Mammograms Alone and in Combination with Tomosynthesis Images." *Radiology*. 2014 Jun;271(3):664-71. Epub 2014 Jan 21. **4** Bernardi D, Macaskill P, Pellegrini M, et al. "Breast cancer screening with tomosynthesis (3D mammography) with acquired or synthetic 2D mammography compared with 2D mammography alone (STORM-2): a population-based prospective study." *Lancet Onc*. 2016 Aug;17(8):1105-1113. Epub 2016 June 23 **5** Durand M, Raghu M, Geisel J, et al. "Synthesized 2D Mammography + Tomosynthesis: Can We See Clearly?" (paper presented at the annual meeting of the Radiological Society of North America, Chicago, IL, December 2015). **6** Choi J, Han B, Ko E, et al. "Comparison with Two-Dimensional Synthetic Mammography Reconstructed from Digital Breast Tomosynthesis and Full Field Digital Mammography for the Detection of T1 Breast Cancer." *European Radiology*. 2016 Aug;26(8):2538-46. Epub 2015 Dec. **7** Woo O, Choi G, Shin H, et al. "Comparative Diagnostic Value of Two-dimensional Synthesized Mammogram and Conventional Full-field Digital Mammogram for Evaluation of Breast Cancer" (poster presented at the annual meeting of the Radiological Society of North America, Chicago, IL, December 2015). **8** Zuckerman S, Conant E, Keller B, et al. "Implementation of Synthesized Two-dimensional Mammography in a Population-based Digital Breast Tomosynthesis Screening Program." *Radiology*. 2016 Dec;281(3):730-736. Epub 2016 July 28. **9** <http://www.fda.gov/downloads/AdvisoryCommittees/CommitteesMeetingMaterials/MedicalDevices/MedicalDevicesAdvisoryCommittee/RadiologicalDevicesPanel/>

SS-00533-EUR-EN Rev.001 (7/17) Hologic Inc. ©2017 All rights reserved. Hologic, 3D, 3Dimensions, 3D Mammography, C-View, Dimensions, Selenia, The Science of Sure, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

[www.hologic.com](http://www.hologic.com) | [info@hologic.com](mailto:info@hologic.com) | + 32.2.711.4680

**No Compromise. No Comparison.**