

HOLOGIC

Better Together™

MyoSure® + Fluent® Pro System

We believe Better is Possible™

Hologic continues to prioritize women's health through partnership and innovation of technologies in the operating room and outpatient suite. Your hysteroscopic solutions now benefit from our improved fluid management system: **Fluent® Pro with SmartFlo™ technology.**

Making better possible, together

The Fluent Pro system was designed to optimize your MyoSure hysteroscopy procedures by increasing efficiency, confidence and performance with **better suction, visualization and resection.**



Hysteroscopy and Fluid Management: MyoSure® + Fluent® Pro System




BETTER SUCTION

Up to 4x faster suction*

Achieve 4x faster suction using the high suction setting on the Fluent Pro system.*

Select between three different levels of suction for resection of simple-to-complex pathology and notice faster draining from the under buttocks drape for end-of-procedure results in less time.




BETTER VISUALIZATION†

Powered by SmartFlo technology

Optimize visualization with SmartFlo technology, featuring adaptive learning and proactive inflow response.

This innovative technology collects inflow rate data and proactively adjusts to deliver consistent pressure during each procedure.



BETTER RESECTION

Up to 2.5x more tissue

Resect 2.5x more tissue with the MyoSure XL device, 2x more tissue with the MyoSure Lite device and 1.7x more tissue with the MyoSure Reach device using the high suction setting.*

The Fluent Pro system implements a 10mmHg pressure setpoint decrease during cutting. This, combined with adjustable suction levels, allows for more tissue to be pulled into the cutting window.

Better Together

Experience the power of the
MyoSure + Fluent Pro system

*Based on a comparison of bench testing performed on the Fluent Pro system high suction setting and the Fluent system default/non-adjustable setting, DHM-16541.

†Based on a survey of 34 physician users during Fluent Pro Limited Market Release, 2024, DHM-16542.

Ordering information		
Product description	Product number	Quantity
TISSUE REMOVAL		
MyoSure Lite Tissue Removal Device	30-403LITE	3-pack
MyoSure Reach Tissue Removal Device	10-403FC	3-pack
MyoSure XL Tissue Removal Device	50-603XL	3-pack
FLUID MANAGEMENT		
Fluent Pro Fluid Management System	FLT-200	1 system
Fluent Pro Disposable Packs	FLT-212	6-pack

Contact your Hologic sales representative for more information.

SS-01636-CAN-001 Rev 001 © 2024 Hologic Inc. All rights reserved. Hologic, MyoSure®, Fluent® Pro, Omni®, and The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.